

EVA Series Approval Label



Bullard
 1898 Safety Way
 CYNTHIANA, KY 41031 USA
 877-BULLARD (285-5273)



MODEL EVA/EVAHL SERIES POWERED AIR PURIFYING RESPIRATORS
 This respirator is approved only in the following configurations:

RESPIRATOR COMPONENTS		PROTECTION1	ALTERNATE HOOD ASSEMBLIES	ALTERNATE SUSPENSION / HARD HAT ASSEMBLIES	ALTERNATE BREATHING TUBES	ALTERNATE BLOWER UNIT ASSEMBLIES	FILTER ASSEMBLY	ALTERNATE CARTRIDGE ASSEMBLY	ALTERNATE BATTERY CONFIGURATIONS	ALTERNATE BELT ASSEMBLIES	ALTERNATE CHARGER ASSEMBLIES	ACCESSORIES	CAUTIONS AND LIMITATIONS2
21C-0836	HE		20LFXL										
21C-0972	HE		20LFL										
23C-2510	OV/GD/CL/HC/HF/SD/HE		20LFXL										
23C-2765	OV/GD/CL/HC/HF/SD/HE		20LFL										
23C-2886	AM/CL/HC/MA/FM/HF/SD/CD/HE		20LFXL										
23C-3097	AM/CL/HC/MA/FM/HF/SD/CD/HE		20LFL										

1 Protection

HE - High Efficiency Particulate Air Filter for Powered, Air Purifying Respirators

AM - Ammonia

CD - Chlorine dioxide

CL - Chlorine

FM - Formaldehyde

HC - Hydrogen chloride

HF - Hydrogen fluoride

MA - Methylamine

SD - Sulfur Dioxide

OV - Organic Vapor

2 Cautions and Limitations

A. Not for use in atmospheres containing less than 19.5% oxygen.

B. Not for use in atmospheres immediately dangerous to life or health.

C. Do not exceed maximum use concentrations established by regulatory standards.

F. Do not use this respirator if airflow is less than four cfm (115 lpm) for tight-fitting facepieces or six cfm (170 lpm) for hoods and / or helmets.

H. Follow established cartridge and canister change schedules or observe ESILI to ensure that cartridges and canisters are replaced before breakthrough occurs.

I. Contains electrical parts that may cause an ignition in flammable or explosive atmospheres.

J. Failure to properly use and maintain this product could result in injury or death.

L. Follow the manufacturer's instructions for changing cartridges and / or filters.

M. All approved respirators shall be selected, fitted, used and maintained in accordance with MSHA, OSHA and other applicable regulations.

N. Never substitute, modify, add or omit parts. Use only exact replacement parts in the configuration specified by the manufacturer.

O. Refer to User's Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.

P. NIOSH does not evaluate respirators for use as surgical masks.